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English

Type

Electric micromotor hose with revolving connector for use with the MC2, the Isolite 100 and 300, with or without light. Quality hoses straight or coiled, in black, grey or white. Standard length 1.7 m or special lengths up to 3 m.

Intended use

Product intended for professional use only. Use in dentistry for prophylaxis, general dentistry and endodontic work.

Technical data and assembly

Classification

Class IIa, «active medical devices», as per section III «Classification», Rule 2 of Annex IX to European Directive 93/42/EEC, namely: «Non-invasive devices intended for channelling liquids or gases for the purpose of administration into the body, which may be connected to an active medical device in Class IIa».

Fitting

Keep the initial alignment of the leads and conduits unchanged.

Place the back-up ring **8** in the locking area **fig. 1**.

Tie the security string **9 fig. 2**.

Description fig. 3 - 4

- 1 Ø 1.5/2.5 mm air spray (small white) with "A" or (blue)
- 2 Ø 1.5/2.5 mm water spray (small white) with "W" or (green)
- 3 Ø 1.5/2.5 mm motor cooling (small white) or (transparent) Ø 2.8/4.1 mm
- 4 (+) red: motor
- 5 black: motor
- 6 (+) brown: bulb
- 7 (0V) blue: bulb
- 8 back-up ring
- 9 security string (if existing)

Accessories

O-Ring joint REF 1300414-010. Changing the O-ring joint **fig. 5**. The split ring REF 212.65.01 **fig. 6** is required only for connection to an

MC2 IR or GTAV motor **fig. 7-8**.

Connection (without the split ring) to an Isolite electric motor **fig. 9-10**. The conical piece **fig. 11** can be useful for connecting the cooling tube to the Unit tube **7 fig. 12**.

Maintenance

Non-sterilisable

Cleaning

Carefully clean all surfaces with a cloth impregnated with Spraynet, **fig. 13**.

Disinfection

Disinfect using a clean cloth soaked in a suitable product. Products containing acetone, chlorine and bleaches are not recommended as disinfectants. To keep the surface of the hose in good condition, it is advisable to periodically wipe the complete length of it with a cloth dusted with talcum powder.

Never immerse in disinfectant solutions.

Not suitable for ultrasonic bath.

Transport

Temperature between -40°C (-40°F) and 70°C (158°F), relative humidity between 10% and 100%, atmospheric pressure 50 kPa to 106 kPa (7.3 to 15.3 psi).

Information

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact Bien-Air Dental SA at the address indicated on the back cover.

Other precautions for use

The device must be used by a qualified person in accordance with the current legal provisions concerning industrial safety, health and accident prevention measures, and these working instructions. In accordance with these requirements, the operators:

- must only use operating devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by Bien-Air Dental;
- must ensure that the device is used only for the purpose for which it is intended, must protect themselves, their patients and third parties from any danger, and must avoid contamination through the use of the product.

The device is intended for medical treatment only; any use other than

that for which this product is intended is unauthorised and may be dangerous. The medical device meets all the current legal requirements.

The device is not authorised for use in an explosive atmosphere (anaesthetic gas).

Recommendations

It is essential to use dry, purified compressed air in order to ensure the long working life of the device. Maintain the quality of the air and the water by regular maintenance of the compressor and the filtration systems. The use of unfiltered hard water will lead to early blockage of the tubes, connectors and spray cones.

Rest the device on a suitable support to avoid risks of infection for yourself, the patient or third parties.

Only use maintenance products and components from Bien-Air Dental. The use of other products and components can void the guarantee.

Guarantee

Terms of guarantee

Bien-Air Dental grants the user a guarantee covering all functional defects, material or production faults. The device is covered by this guarantee for 12 months from the date of invoicing.

Bien-Air Dental offers a 24-month guarantee for the glass-bar light conductors.

In the event of justified claim, Bien-Air Dental or its authorised representative will fulfil the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

Bien-Air Dental shall not be held responsible for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear
- improper use

- non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences
- poor connections, whether of the air, water or electricity supply.

The guarantee does not cover flexible "fibre optic" type conductors, or any parts made of synthetic materials.

The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorised by Bien-Air Dental.

Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

Designation micromotor		MC2 40 IR	MC2 40 GTAV	MC2 Isolite 300	MC2 Isolite 100
REF micromotor		1600073-001	1600075-001	1600078-001	1600079-001
Designation hose	REF hose				
Grey hose with fixed connector	1600120-001	x	(x)	x	x
Black hose with fixed connector	1600297-001	x	(x)	x	x
Grey hose, revolving connector	1600298-001			x	x
Black hose, revolving connector	1600299-001			x	x
White hose with fixed connector	1600314-001	x	(x)	x	x
Black coiled hose with fixed connector	1600315-001	x	(x)	x	x
White hose, revolving connector	1600316-001			x	x

(X) Do not connect supply of water and air spray (risk to damage the micromotor)

REF	Legend	List of Bien-Air Dental SA registered trade mark products ®:				
1600120-001	Grey hose with fixed connector, for micromotors MC2 and Isolite	Aquilon®	ChiroPro®	Gyrolina®	Lubrimed®	Prestilina®
1600297-001	Black hose with fixed connector, for micromotors MC2 and Isolite	Bora®	Eolia®	Isolite®	MX®	Spraynet®
1600298-001	Grey hose, revolving connector for micromotors Isolite	Boralina®	Gyro®	Lubrifiuid®	PowerCare®	
1600299-001	Black hose, revolving connector for micromotors Isolite	In these instructions, “Device” corresponds to the product described in the heading “Type”. For example, turbine, contra-angle, handpiece, micromotor, tube, electronics, connectors, station etc.				
1600314-001	White hose with fixed connector, for micromotors MC2 and Isolite					
1600315-001	Black coiled hose with fixed connector for micromotors MC2 and Isolite					
1600316-001	White hose, revolving connector for micromotors Isolite	Symbols				
1300414-010	O-Ring joint	 Manufacturer.	 CE Marking with number of the notified body.			
212.65.01-001	Split ring	 Wear rubber gloves.	 Move in the direction indicated.			
1600036-006	Spraynet, cleaning spray 500 ml, box of 6 cans	 Attention.	 Move fully to the stop, in the direction indicated.			

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